

CLONING MIST

KEEP OUT OF REACH OF CHILDREN SHAKE WELL BEFORE USING

Please see Rx Green Technologies Growing Schedules for more specific application rates

WARRANTY: Rx Green Technologies, LLC warrants that this product conforms to the analysis on its label. When used in accordance with label directions, under normal conditions, this product is reasonably fit for its intended purposes. Since timing, method of application, weather, plant and soil conditions, mixture with other chemicals and other factors affecting the use of this product are beyond our control, no warranty is given to the use of this product contrary to label directions or under conditions which are abnormal or not reasonably foreseeable. The user assumes all risks of any such case.

LIFE™

CLONING MIST

RX GREEN
TECHNOLOGIES

12 FLUID OUNCES (355ML)
NET WEIGHT: 0.78 LBS (0.35 KG)

CLONING MIST

LIFE™ Cloning Mist is an innovative ready-to-use spray for cuttings.



FOR BEST RESULTS USE WITH:

- LIFE™ Cloning Solution
- LIFE™ Cloning Gel

DIRECTIONS FOR USE:

STEP 1: Spray the leaves of stock plants with LIFE Cloning Mist one week prior to taking cuttings.

STEP 2: Take cuttings from stock plant and place into growing medium. (aeroponics/rockwool/peat/coco)

STEP 3: Spray cuttings after planted in growing medium.

STEP 4: Spray cuttings every 2-3 days.

APPLY ONLY AS DIRECTED

RX GREEN
TECHNOLOGIES

Rx Green Technologies, LLC
www.rxgreentechnologies.com
PO Box 3213 Manchester, NH 03105



Lot Code: